

M1050F3



### Description

This Thorlabs Fiber-Coupled LED is mounted to the end of a heat sink. The output is compatible with SMA fiber connectors. This LED needs to be supplied with a constant current that must not exceed the maximum current. The current source must be able to deliver this current at the specified forward voltage. For specifications, please see the table below.

### Specifications

Specification	Value
Color	IR
Nominal Wavelength	1050 nm
Bandwidth (FWHM)	37 nm
Emitter Size	1 mm x 1 mm
Maximum Current (CW)	600 mA
Electrical Power	840 mW
Typical Lifetime	>10 000 h
Operating Temperature (Non-Condensing)	0 to 40 °C
Storage Temperature	-40 to 70 °C
Risk Group <sup>a</sup>	RG0 - Exempt Group

a. According to the standard IEC 62471:2006, Photobiological Safety of Lamps and Lamp Systems.

	Symbol	Min	Typical	Max
Peak Wavelength <sup>a</sup>	$\lambda_p$	1000 nm	1050 nm	1100 nm
LED Output Power (Ø400 µm Fiber) <sup>b,c</sup>	$P_{out}$	2.3 mW	3.0 mW	-
LED Output Power (Ø200 µm Fiber) <sup>b,d</sup>	$P_{out}$	-	0.92 mW	-
Forward Voltage <sup>b</sup>	$V_F$	-	1.4 V	1.6 V

a. When Driven with Maximum Current; Valid for Bare Die

b. When Driven with the Maximum Current

c. For multimode fiber with a Ø400 µm core and 0.39 NA (Item # FT400EMT).

d. For multimode fiber with a Ø200 µm core and 0.22 NA (Item # FG200UCC).

### Operating Instructions

Be sure to provide air ventilation in order to avoid overheating, drops in optical power, and reduced lifetime. Be aware that each LED has a characteristic switch-on behavior, which depends on the LED properties and environment conditions. To supply sufficient heat dissipation, this LED has a unique thermal design that reduces the power decay to a minimum.

The drawing on page two shows the LED's male connector which is a standard M8x1 sensor circular connector. Pins 1 and 2 are connected to the LED. Pins 3 and 4 are used for the internal EEPROM. This pin assignment is valid for all Thorlabs LED drivers. For use of third party LED drivers, please ensure the correct pin assignment.



## ***Power Supply***

---

Please see the Thorlabs website for the available LED Drivers. When selecting the LED Driver, please ensure that:

- The operating current never exceeds the Maximum Current specified in the LED Specification table.
- Sufficient forward voltage is supplied.

If you decide to use your own DC source, please ensure that the correct connection is made to Pins 1 and 2.

## ***Maintenance and Service***

---

Do not stick any items into the SMA connector aperture - you may damage the LED.

This LED is not water resistant and must be protected from adverse weather conditions. To avoid damage, do not expose it to spray, liquids, or solvents. This LED does not contain any parts serviceable by the user and does not require regular user maintenance. Do not open the enclosure. If a malfunction occurs or if your fiber-coupled LED needs repair, contact Thorlabs for return instructions.

## ***Warnings and Safety***

---

During normal operation, the casing temperature may exceed ambient temperature by as much as 25 °C (45 °F). To prevent higher case temperatures, the product should be operated without anything hindering air movement around the convective cooling fins.

Please note that this product is not suitable for household room illumination.

This LED must not be operated in explosive environments and should only be used with shielded connection cables.

All statements regarding safety of operation and technical data only apply when the unit is operated correctly according to its specifications. The safety of any system incorporating the equipment is the responsibility of the assembler of the system.

Inappropriate use of any Mounted LED product may result in permanent eye damage. To prevent injury, use this product in accordance with the International Standard "Photobiological Safety of Lamps & Lamp Systems" IEC 62471. This LED falls under the Risk Group RG0 - Exempt Group in accordance to the standard IEC 62471:2006.

If using this LED in a microscope application as a replacement for mercury vapor lamp, the same precautions should be taken.

### **IR Warning Statement**

This LED emits intense IR radiation during operation. Precautions must be taken to prevent looking directly at the IR light. If viewing the IR light is necessary, IR light protective glasses must be worn to avoid eye damage. Do not look directly into the IR light or look through the optical system during operation, as this can be harmful to the eyes, even for brief periods, due to the high intensity of the IR light.

